

# Future Prospects with DKIST

Tom Schad – DLNIRSP Instrument Scientist – Univ. of Hawaii/IfA ISSI Workshop on Coronal Rain: 23 – 27 February 2015, Bern, Switzerland







# Senator Daniel K. Inouye Solar Telescope (DKIST) [previously ATST]



- 4 m aperture
- Haleakala, HI (10010 ft; 3050 m)
- Gregorian, off-axis, Alt-Az
- Coude instruments:
  - VBI
  - VISP
  - ► VTF
  - DL-NIRSP
  - CRYO-NIRSP
- FOV: 2.8 (max 5) arcmin

- DKIST will be operational and ready to carry out its Critical Science Plan (CSP) in Fall 2019.
- What critical coronal rain observations can our team propose for DKIST?
  - Target for next meeting: draft experiment details using designed instrument capabilities?



dkist.nso.edu/science

## Three reasons we need DKIST for coronal rain

### #I Coronal structure is not resolved!



## Three reasons we need DKIST for coronal rain

- #I Coronal structure is not resolved!
- #2 Photons!
  - Rain signals very weak in optical/IR diagnostics
  - Where is the rest of the coronal rain mass?
  - Time-resolved chromospheric polarimetry
- #3 Coronal Field Diagnostics (esp. in near/mid IR)
  - Control of scattered light  $\rightarrow$  Off-axis design
  - Active occulter at Gregorian focus; Lyot stop
  - No windows! -> Extends max wavelength range out to ~ 28 um

# **DKIST Construction Update**

- Construction started on Haleakala late 2012
- Mirror fab started well in advance of 2012.
- Large systems (i.e. telescope mount, enclosure, etc.) near/ past acceptance in many cases.
- Now the presence on the mountain is very real!

### Stress-lap polishing of M1 commissioning blank (U. of Arizona)



Movie courtesy of NSO/DKIST

Telescope enclosure assembled in Spain by IDOM.

Will be installed on site this spring!



### Telescope Mount Assembly

Image courtesy NSO/DKIST

### Telescope Mount Assembly – Slewing



Movie courtesy of NSO/DKIST



## Coude Lab Platform and Mount in Factory



Movie courtesy of NSO/DKIST and Ingersoll

16 February 2015 Photo by Tom Schad STO TOW

FFF

## **DKIST Facility Instruments**





## DKIST Spectral Distribution (i.e. the Autobahn)

- Spectral distribution controlled by series of dichroic beam-splitters downstream of AO DM.
- Manual setup; must be selected for given experiment
- Cryo-NIRSP operates alone; No AO.



#### And now, a **very** brief introduction to the instruments...

VBI VISP VTF CRYO-NISRP DL-NIRS	P
---------------------------------	---

# VBI: The Visible Broadband Imager

- Developed by the National Solar Observatory (NSO), Boulder CO
- Fast 4K x 4K SPECKLE-imaging! Up to 3.2 sec cadence for reconstructed images
   VBI H-Beta Filter Obs at DST
- Blue-red channels simultaneously; synchronized to 10 ms
- FOV: Up to 2' square;
  45" / 69" physical (blue/red)
- II and I7 mas sampling!
- VBI-blue:
  - <u>393 nm (rain!)</u>
  - <u>H-beta 430 nm (rain!)</u>
  - 450 nm
  - ▶ 486 nm
- VBI-red:
  - <u>656.3 nm (rain!)</u>
  - ▶ 668.4 nm
  - ▶ 705 nm
  - ? (hot coronal filter?)



#### Movie courtesy of NSO/DKIST

Y₿I

NIRSF

**DL-NIRSP** 

# VISP: The Visible Spectropolarimeter

ViSP OVERVIEW

- Developed by High Altitude Observatory (HAO), Boulder, CO
- ▶ 380 to 900 nm [380 1600 nm goal]
- Up to three line simultaneously
- ▶ 0.07" spatial resolution
- ▶ R >= 180000
- I0^-3 I<sub>cont</sub> polarimetry in 10 sec

An extraordinarily flexible slit-spectropolarimeter!

For coronal rain: Sit-and-stare polarimetry?

VTF

**NIRSF** 

**VBI** 

# VTF: Visible Tunable Filter

- Developed by KIS, Freiburg, Germany
- Fabry-perot imaging spectropolarimeter; similar in functionality to CRISP!
- Rain: He D3; H alpha; Ca II 854 nm

Description	Value	Comment
Wavelength range	520 – 870 nm	
Spatial resolution	0.028"	Diffraction limited
Spectral Resolution	6 pm	at 600 nm
Spectral sampling	3pm	at 600 nm
Cadence, Full Stokes	13 s	SNR=650, 30 Hz frame rate
Field of view	60" x 60 "	4 k x 4k detector
Doppler sensitivity	85 m/s	
Polarimetric sensitivity	3 x 10 <sup>-3</sup>	SNR=650
Magnetic field sensitivity	20 G	longitudinal



Y₿I

VISP

# Cryo-NIRSP: The Cryogenic Near-IR Spectropolarimeter

- Developed by University of Hawaii Institute for Astronomy (IfA)
- Operates alone; No AO
- Huge FOV (5'); deep cyro-genic instrument. Specializing in coronal magnetic fields!
- > 500 to 5000 nm
- First set of lines:
  - Fe XIV 530 nm
  - Fe XIII 1074.7
  - He I 1083 nm
  - ▶ Si X 1430 nm
  - ▶ Si IX 3930 nm
  - CO 4651 nm



# The Diffraction-Limited Near IR Spectropolarimeter (DL-NIRSP)

- Developed by University of Hawaii Institute for Astronomy (IfA)
- The first IFU based spectropolarimeter for solar physics!



NISPR

**DL-NIRSP** 

## The Diffraction-Limited Near IR Spectropolarimeter (DL-NIRSP)



One tile at a time, DL-NIRSP builds spectropolarimetric full data cubes: [X;Y; $\lambda$ ; S [=I,Q,U,V];t]

Y₿I

CRYO.

VISP

VTF

CRYO-NISPR

## **Quick Facts**

Instrument type:	Reflection-grating based, multi-arm integral field spectrograph			
Spectral Range:	500 nm to 1800 nm			
Number of cameras/channels:	3			
Spectral channel ranges:	500 nm to 900 nm (most subject to beam sharing!)			
	900 nm to 1350 nm			
	1350 nm to 1800 nm			
Field of View:	Up to 165" circular (120" x 120" square)			
Polarimetry:	Full Stokes (or Spectroscopy-only)			
Target Polarimetric accuracy:	> 5 x 10 <sup>-4</sup> I <sub>c</sub> (dependent on polarization amplitudes)			
Parallel Operations:	Can observed in parallel with WFC, VBI, VISP, VTP			
Target domain:	On-disk and off-limb			

## DL-NIRSP's enabling technology

BiFOIS: The Birefringent Fiber-Optic Image Slicer



Lin & Versteegh SPIE (2006); Schad et al. SPIE (2014)

## Real examples of prototype spectra

·								- Similarity		-		-
	COMPANY OF THE OWNER.									STATISTICS.	-	-
				망가 문제			Contraction of the		-	PROVIDE NO.		
and the second sec								-		-	-	
							10000	-	-	and the second		-
the second se			and the second sec	나라 것 같아.			ALC: NO.	-	And in case of the local division of the loc	Internet -		-
	THE REAL PROPERTY AND INCOME.							The second secon	and the second second	101000000		and the second second
			CONTRACTOR OF CONTRACTOR				The second second		COLUMN ST	The second second		
											Contraction of the local division of the loc	The second second
			Contraction of the local division of the loc		-			STATISTICS IN				-
											Statistics of the	
				1998년 1998년 <b>-</b>			-				and the second s	-
1000 COL 001											-	
1000	100				. 11			100.5200/100		STATISTICS.		
100 million (100 m		TAXABLE INC.							100 million and	121201		100
						A		States and and			The second second	
	and the second se								and the second		and the second	-
								and in the second	and the second s			and the second se
COMPANY OF A								-		220401025		
			and the second se					the second second			-	
CONTRACTOR DOCUMENTS	Contract Name								200322010	COLUMN IN	Transmission of	The subscription of the su
the second se		100 T				-	-	-				-
Transfer and the second									100818		States and states of	-
							-	100			Statement of the local division of the local	
					-		1.			10/00/112		Contraction of the local division of the loc
1000				- 1 - 1 - 1 - <del>-</del>								And in case of the local division of the loc
1000	and the second second									Sector States and States	and the second second	And in case of the local division of the loc
	The second						Section and an	State of the local division of the local div			COLUMN TWO IS NOT	
						1000	STATISTICS.	100 100 100	-	-	The second second	Concession of the
CONTRACTOR OF CONTRACTOR	and the second second				-		-	No. of Concession, Name	1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.1.	THE REAL PROPERTY OF		
											-	
Contraction of the local division of the loc										Concession of the local division of the loca		-
									and the second s	COLUMN TO BE		-
					initiation in		Control Street Street Street				COLUMN DO UND	
							-	-			THE OWNER WHEN	and the second s
							1.1.1.1.1					
			and the second sec				-	The second se			A DESCRIPTION OF	
	1000							States of the local division of the local di				
							- NY ING AND					

He I 1083 nm

D

Fe I 1565 nm

# SPIES He I 1083 nm





Figure 1: The optical system of the DL-NIRSP, including DKIST M8, M9, and M10



### What will DL-NIRSP spectra look like?

Þ

Synthetic DL-NIRSP Spectral Image - 1082.9 nm



# **DL-NIRSP Spatial Resolution Modes**

Table 1: DL-NIRSP Spatial and Spectral Sampling Resolution Modes (reproduction of DL-NIRSP DRD Table 1).

Mode	Feed Optics	IFU	Spatial Sampling <sup>1</sup>	IFU Format	Unit IFU FOV	Spectral Sampling
High-Res	F/62	BiFOIS- 5K-36	0.03"	80 x 60	2. 4" x 1.8"	~ 250,000
Mid-Res	F/24	BiFOIS- 5K-36	0.077"	80 x 60	6.16" x 4.62"	~ 250,000
Wide-Field	F/24	BiFOIS- 2.5K-72	0.464″	60 x 40	27.84" x 18.56"	~ 250,000

### Summary

- DKIST will be an ordinarily flexible instrument!
- Fast high resolution instrumentation!
- Deep spectropolarimetry of fine scaled structure!







## A shopping list of DKIST rain observations...

- VBI-blue: Ca H/K + H-Beta + Green line?
- VISP: (not initially used for rain; or sit-n-state mode)
- VTF: H-alpha + Ca 854.2 nm
- DL-NIRSP:
  - He I 1083 nm
  - Si X 1430 nm
  - Fe XIII 1074.7 nm
  - Fe XIII 1080 nm

- Cryo-NIRSP
  - **TBD..**